# Comparison of Bioinformatics Pipelines to Analyze Fish Methylomes in Relation to Environmental Adaptation



Mirella Arguelles<sup>1</sup>, Keem Uarren Eiryll B. Antonio<sup>2\*</sup>, Mylah Rystie U. Anacleto<sup>1</sup>, and Jorge Gil C. Angeles<sup>2</sup>

<sup>1</sup>Institute of Computer Science, College of Arts and Sciences <sup>2</sup>Philippine Genome Center - Program for Agriculture, Livestock, Forestry, and Fisheries, University of the Philippines Los Baños \*Corresponding author: kbantonio1@up.edu.ph

## **ABSTRACT**


Methylation of DNA is a well-studied epigenetic mechanism and plays an important role in gene regulation of organisms. Through the recent advancements of next-generation sequencing (NGS) technologies, generation of high-throughput methylation data, and elucidation of DNA methylation affects the organism's adaptation, such as with fishes, to its environment could be attained. Recognition of methylation patterns through NGS and the subsequent detection of significant methylation differences are some challenges to effectively analyze fish methylomes. There are several bioinformatics pipelines employed to analyze methylomes in mammals, insects, and plants but none have been comparatively used to analyze fish methylomes especially to understand their environmental adaptation. In this study, three commonly-used methylation bioinformatic pipelines (Bismark, BiSulfite Bolt, and BWA-Meth) were tested on the sole available methylome data relating to the environmental adaptation of fishes, that from the three-spined stickleback (Gasterosteus aculeatus). Comparative analyses focusing on the performance and functionality of these three programs were conducted to detect each program's advantages and limitations. Bismark and BiSulfite Bolt have no significant differences relating to their run time and are significantly faster than BWA-Meth to analyze the given fish methylome data. Comparing the primary features of the three programs, Bismark has the further advantage to remove duplicated reads given different biological samples. It also has features that would greatly improve data visualization and provide summary reports that can be used for downstream analyses. The findings of this study could be adapted for future methylome analyses concerning the environmental adaptation of fishes.

#### INTRODUCTION


Epigenetics is the study of changes in gene expression that do not involve alterations to the DNA sequence. DNA methylation involves the addition of a methyl group at the 5th carbon of cytosine. It is the most widely studied epigenetic mark and affects a variety of biological processes of an organism such as its development and adaptation to its environment. Through next-generation sequencing techniques, bisulfite modification is employed to study changes of DNA methylation patterns. Several bioinformatics pipelines which specialize in understanding bisulfite sequencing data have been used to identify such methylation patterns. These pipelines have been used in various literature to determine which is most suited for mammals, insects, and plants. However, none have been studied for fish methylomes especially to understand their environmental adaptation.

This study conducts a comparative analysis on commonly used methylome analysis pipelines using available methylomic reads from the three-spined stickleback (Gasterosteus aculeatus). The comparison between the tools were based on run time and evaluation of their primary features. Identifying a suitable pipeline for the project may improve the accuracy and efficiency of DNA methylation analysis for fish methylomes. The results provide researchers with adequate direction on which methylation analysis pipeline to choose for their studies.

# **METHODOLOGY**



### **RESULTS**



- Figure 1. Normalized run times (n = 3) of three different methylation analysis
- No significant differences between Bismark and BiSulfite Bolt ( $\alpha$  < 0.01)
- Both Bismark and BiSulfite Bolt have **faster run time** than BWA-Meth ( $\alpha$  < 0.01)
- Table 1. Summary of key features of three different methylation analysis pipelines


letter have no statistical difference ( $\alpha$  < 0.01).

pipelines using available methylome data relating to the environmental adaptation

of three-spined stickleback (Gasterosteus aculeatus). Software with the same

Feature	Bismark	BiSulfite Bolt	BWA-Meth and MethylDackel
Reference genome indexing	✓	✓	✓
Alignment of bisulfite sequencing reads	✓	✓	✓
Deduplication of sequencing reads	✓	X	X
Methylation calling	✓	✓	✓
Processing report generation	✓	X	✓
Summary report generation	✓	X	X
HTML report visualization	✓	X	X
Python API availability	X	✓	X
Standalone software	✓	✓	X

- Bismark has features for deduplication of sequencing reads and plot visualization of reports while these features do not exist in BiSulfite Bolt and BWA-Meth.
- BWA-Meth is not a standalone methylation analysis software as it requires MethylDackel to perform methylation calls and generate analysis reports.



	track type='	'bedGraph"	descript	tion="SR	R166766	40.sorted	CpG	methylation	levels"
	NC_053212.1	4078	4079	100	1	0			
	NC_053212.1	4087	4088	100	1	0			
	NC_053212.1	4099	4100	100	1	0			
port	NC_053212.1	4104	4105	100	1	0			
3 2022	NC_053212.1	4111	4112	100	1	0			
	NC_053212.1	4115	4116	100	1	0			
	NC_053212.1	4120	4121	100	1	0			
	NC_053212.1	4123	4124	100	1	0			
	NC_053212.1	4135	4136	100	1	0			
	NC_053212.1	4148	4149	100	1	0			
	NC_053212.1	4170	4171	100	1	0			
	NC_053212.1	6321	6322	100	1	0			
	NC_053212.1	6358	6359	100	1	0			
	NC_053212.1	7500	7501	100	1	0			
	NC_053212.1	7549	7550	100	2	0			
	NC_053212.1	7557	7558	100	2	0			
	NC_053212.1	7562	7563	100	2	0			
	NC_053212.1	7683	7684	100	2	0			
	NC_053212.1	7701	7702	100	2	0			
	NC_053212.1	7713	7714	100	2	0			
	NC_053212.1	7717	7718	100	2	0			(b)
	NC_053212.1	7886	7887	100	1	0			

Figure 2. Methylation report outputs of (a) Bismark and (b) BiSulfite Bolt

## CONCLUSION

From the benchmarking tests performed using fish methylome data, it was found out that Bismark and BiSulfite Bolt had the shortest time performance and both were faster than BWA-Meth. Bismark had more comprehensive features in comparison to BiSulfite Bolt and BWA-Meth. The results show that Bismark is a suitable candidate to analyze fish methylomes and other methylomes relating to environmental adaptation.

## SELECTED REFERENCES

